

Tiragolumab ELISA Kit

SKU: GXEL00229 | Assay: Quantitative ELISA (Colorimetric)

The Tiragolumab ELISA Kit (SKU: GXEL00229) is a high-sensitivity quantitative ELISA for the quantitative measurement of Tiragolumab in Plasma, Serum. The assay covers a detection range of 31.25 ng/mL - 2000 ng/mL and is read colorimetrically at 450 nm.

Background: Tiragolumab is an investigational, fully human monoclonal IgG1 antibody designed to prevent PVR-TIGIT interaction.

Specifications

Product Name	Tiragolumab ELISA Kit
SKU	GXEL00229
Accession	Q495A1
Alternative Names	MTIG-7192-A, MTIG-7192A, MTIG7192A, RG-6058, RO-7092284, RO7092284, 1918185-84-8
Assay Type	Quantitative
Detection Method	Colorimetric
Detection Range	31.25 ng/mL - 2000 ng/mL
Recovery	75-110%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at -20°C, the rest reagents should be store at 4°C. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00229.pdf> | Version 1.0, 09 July 2026