

Frovocimab ELISA Kit

SKU: GXEL00234 | Assay: Quantitative ELISA (Colorimetric)

The Frovocimab ELISA Kit (SKU: GXEL00234) is a high-sensitivity quantitative ELISA for the quantitative measurement of Frovocimab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

Background: Frovocimab (LY3015014) (INN) is a monoclonal antibody designed for the treatment of hypercholesterolemia. This drug is being developed by Eli Lilly and Company.

Specifications

Product Name	Frovocimab ELISA Kit
SKU	GXEL00234
Accession	Q8NBP7
Alternative Names	LY-3015014, 1643672-70-1
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	0.156 µg/ml
Detection Range	0.31-5 µg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00234.pdf> | Version 1.0, 09 July 2026