

## Tezepelumab ELISA Kit

SKU: **GXEL00238** | Assay: Quantitative ELISA (Colorimetric)

The Tezepelumab ELISA Kit (SKU: GXEL00238) is a high-sensitivity quantitative ELISA for the quantitative measurement of Tezepelumab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

**Background:** Tezepelumab is a first-in-class human monoclonal antibody that binds to TSLP, thus inhibiting its interaction with TSLP receptor complex. Tezepelumab given as an add-on-therapy to patients with severe uncontrolled asthma has shown safety, tolerability and efficacy. Several trials are evaluating the long-term safety and the efficacy of tezepelumab in adults and adolescents with severe uncontrolled asthma.

### Specifications

<b>Product Name</b>	Tezepelumab ELISA Kit
<b>SKU</b>	GXEL00238
<b>Accession</b>	Q969D9
<b>Alternative Names</b>	AMG 157, AMG-157, MEDI-9929, MEDI9929, anti-TSLP-Receptor (TSLP-R, CrI2), 1572943-04-4
<b>Assay Type</b>	Quantitative
<b>Detection Method</b>	Colorimetric
<b>Sensitivity</b>	0.156 µg/ml
<b>Detection Range</b>	0.31-5 µg/mL
<b>Recovery</b>	80-120%
<b>Sample Types</b>	Plasma, Serum
<b>Storage</b>	2-8°C. See Storage & Stability.
<b>Note</b>	For Research Use Only. Not for diagnostic procedures.

### Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | [www.assaygenie.com](http://www.assaygenie.com) | [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com)  
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00238.pdf> | Version 1.0, 09 July 2026