

Ianalumab ELISA Kit

SKU: GXEL00241 | Assay: Quantitative ELISA (Colorimetric)

The Ianalumab ELISA Kit (SKU: GXEL00241) is a high-sensitivity quantitative ELISA for the quantitative measurement of Ianalumab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

Background: Ianalumab (VAY736) is an anti-B-cell activating factor (BAFF) receptor fully human IgG1 monoclonal antibody, engineered for direct ADCC-mediated B-cell depletion, thus providing a dual mode of action and targeted approach to treat pSS. The current phase 2b study aimed at establishing a dose-response relationship over a wide range of VAY736 doses, using change from baseline (BL) in the EULAR Sjogren's Syndrome Disease Activity index (ESSDAI) over 24 Weeks (Wks) as the primary endpoint. Here we report Wk24 efficacy and safety. The study is ongoing with a second blinded treatment period up to Wk52.

Specifications

Product Name	Ianalumab ELISA Kit
SKU	GXEL00241
Accession	Q96RJ3
Alternative Names	VAY-736, 1929549-92-7
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	0.156 µg/ml
Detection Range	0.31-5 µg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lhx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00241.pdf> | Version 1.0, 09 July 2026