

Lirilumab ELISA Kit

SKU: GXEL00262 | Assay: Quantitative ELISA (Colorimetric)

The Lirilumab ELISA Kit (SKU: GXEL00262) is a high-sensitivity quantitative ELISA for the quantitative measurement of Lirilumab in Plasma, Serum. The assay covers a detection range of 0.31-5 µg/mL with a sensitivity of 0.156 µg/ml and is read colorimetrically at 450 nm.

Background: Lirilumab (IPH2102/BMS-986015/BMS-986015-01) is a fully human monoclonal antibody (mAb) that is designed to block the interaction between KIR2DL-1, -2, -3 inhibitory receptors and their ligands. This drug was developed by Innate Pharma and is licensed to Bristol-Myers Squibb. A phase 2 clinical trial for acute myeloid leukemia (AML) was terminated early (

Specifications

Product Name	Lirilumab ELISA Kit
SKU	GXEL00262
Alternative Names	BMS-986015, IPH2102, 1000676-41-4
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	0.156 µg/ml
Detection Range	0.31-5 µg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
Datasheet: <https://store-h68l9z2lhx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00262.pdf> | Version 1.0, 09 July 2026