

Anti-SARS-CoV-2 RBD (Delta) Human IgG ELISA Kit

SKU: GXEL00275 | Assay: Quantitative ELISA (Colorimetric)

The Anti-SARS-CoV-2 RBD (Delta) Human IgG ELISA Kit (SKU: GXEL00275) is a high-sensitivity quantitative ELISA for the quantitative measurement of Anti-SARS-CoV-2 RBD (Delta) Human IgG in Plasma, Serum.

Specifications

Product Name	Anti-SARS-CoV-2 RBD (Delta) Human IgG ELISA Kit
SKU	GXEL00275
Assay Type	Quantitative
Detection Method	Colorimetric
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00275.pdf> | Version 1.0, 09 July 2026