

## Anti-SARS-CoV-2 Spike Protein (Trimer) Human IgA ELISA Kit

SKU: GXEL00278 | Assay: Quantitative ELISA (Colorimetric)

The Anti-SARS-CoV-2 Spike Protein (Trimer) Human IgA ELISA Kit (SKU: GXEL00278) is a high-sensitivity quantitative ELISA for the quantitative measurement of Anti-SARS-CoV-2 Spike Protein (Trimer) Human IgA in Plasma, Serum.

### Specifications

|                         |                                                            |
|-------------------------|------------------------------------------------------------|
| <b>Product Name</b>     | Anti-SARS-CoV-2 Spike Protein (Trimer) Human IgA ELISA Kit |
| <b>SKU</b>              | GXEL00278                                                  |
| <b>Assay Type</b>       | Quantitative                                               |
| <b>Detection Method</b> | Colorimetric                                               |
| <b>Recovery</b>         | 80-120%                                                    |
| <b>Sample Types</b>     | Plasma, Serum                                              |
| <b>Storage</b>          | 2-8°C. See Storage & Stability.                            |
| <b>Note</b>             | For Research Use Only. Not for diagnostic procedures.      |

### Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | [www.assaygenie.com](http://www.assaygenie.com) | [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com)  
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00278.pdf> | Version 1.0, 09 July 2026