

Anti-SARS-CoV-2 Spike Protein Mouse IgA ELISA Kit

SKU: GXEL00316 | Assay: Quantitative ELISA (Colorimetric)

The Anti-SARS-CoV-2 Spike Protein Mouse IgA ELISA Kit (SKU: GXEL00316) is a high-sensitivity quantitative ELISA for the quantitative measurement of Anti-SARS-CoV-2 Spike Protein Mouse IgA in Plasma, Serum. The assay covers a detection range of 125 - 8, 000 pg/mL with a sensitivity of 47.47 pg/mL and is read colorimetrically at 450 nm.

Specifications

Product Name	Anti-SARS-CoV-2 Spike Protein Mouse IgA ELISA Kit
SKU	GXEL00316
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	47.47 pg/mL
Detection Range	125 - 8, 000 pg/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00316.pdf> | Version 1.0, 09 July 2026