

## Anti-Rabies Virus (RV) Glycoprotein (G) Dog IgG ELISA Kit

SKU: GXEL00340 | Assay: Quantitative ELISA (Colorimetric)

The Anti-Rabies Virus (RV) Glycoprotein (G) Dog IgG ELISA Kit (SKU: GXEL00340) is a high-sensitivity quantitative ELISA for the quantitative measurement of Anti-Rabies Virus (RV) Glycoprotein (G) Dog IgG in Plasma, Serum. The assay covers a detection range of 1.56 - 100 ng/mL with a sensitivity of 1.09 ng/mL and is read colorimetrically at 450 nm.

### Specifications

|                         |                                                           |
|-------------------------|-----------------------------------------------------------|
| <b>Product Name</b>     | Anti-Rabies Virus (RV) Glycoprotein (G) Dog IgG ELISA Kit |
| <b>SKU</b>              | GXEL00340                                                 |
| <b>Assay Type</b>       | Quantitative                                              |
| <b>Detection Method</b> | Colorimetric                                              |
| <b>Sensitivity</b>      | 1.09 ng/mL                                                |
| <b>Detection Range</b>  | 1.56 - 100 ng/mL                                          |
| <b>Recovery</b>         | 80-120%                                                   |
| <b>Sample Types</b>     | Plasma, Serum                                             |
| <b>Storage</b>          | 2-8°C. See Storage & Stability.                           |
| <b>Note</b>             | For Research Use Only. Not for diagnostic procedures.     |

### Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | [www.assaygenie.com](http://www.assaygenie.com) | [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com)  
Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00340.pdf> | Version 1.0, 09 July 2026