

Insulin Aspart ELISA Kit

SKU: GXEL00376 | Assay: Quantitative ELISA (Colorimetric)

The Insulin Aspart ELISA Kit (SKU: GXEL00376) is a high-sensitivity quantitative ELISA for the quantitative measurement of Insulin Aspart in Plasma, Serum. The assay covers a detection range of 78.125 - 5, 000 ng/mL with a sensitivity of 33.12 ng/mL and is read colorimetrically at 450 nm.

Specifications

Product Name	Insulin Aspart ELISA Kit
SKU	GXEL00376
Assay Type	Quantitative
Detection Method	Colorimetric
Sensitivity	33.12 ng/mL
Detection Range	78.125 - 5, 000 ng/mL
Recovery	80-120%
Sample Types	Plasma, Serum
Storage	2-8°C. See Storage & Stability.
Note	For Research Use Only. Not for diagnostic procedures.

Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00376.pdf> | Version 1.0, 09 July 2026