

## SARS-CoV-2 Surrogate Virus Neutralization Test Kit (JN.1)

SKU: GXEL00381 | Assay: Quantitative ELISA (Colorimetric)

The SARS-CoV-2 Surrogate Virus Neutralization Test Kit (JN.1) (SKU: GXEL00381) is a high-sensitivity quantitative ELISA for the quantitative measurement of SARS-CoV-2 Surrogate Virus Neutralization Test Kit (JN.1) in Plasma, Serum. The assay covers a detection range of 31.25 - 2000 ng/mL with a sensitivity of 16.81 ng/mL and is read colorimetrically at 450 nm.

### Specifications

<b>Product Name</b>	SARS-CoV-2 Surrogate Virus Neutralization Test Kit (JN.1)
<b>SKU</b>	GXEL00381
<b>Assay Type</b>	Quantitative
<b>Detection Method</b>	Colorimetric
<b>Sensitivity</b>	16.81 ng/mL
<b>Detection Range</b>	31.25 - 2000 ng/mL
<b>Recovery</b>	80-120%
<b>Sample Types</b>	Plasma, Serum
<b>Storage</b>	2-8°C. See Storage & Stability.
<b>Note</b>	For Research Use Only. Not for diagnostic procedures.

### Storage & Stability

The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% prior to the expiration date under appropriate storage condition. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | [www.assaygenie.com](http://www.assaygenie.com) | [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com)  
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00381.pdf> | Version 1.0, 09 July 2026