

## Anti-Mazdutide ADA ELISA Kit

SKU: GXEL00395 | Assay: Quantitative ELISA (Colorimetric)

The Anti-Mazdutide ADA ELISA Kit (SKU: GXEL00395) is a high-sensitivity quantitative ELISA for the quantitative measurement of Anti-Mazdutide in Plasma, Serum. The assay covers a detection range of 15.63 - 10, 000 ng/mL with a sensitivity of 6.49 ng/mL and is read colorimetrically at 450 nm.

### Specifications

<b>Product Name</b>	Anti-Mazdutide ADA ELISA Kit
<b>SKU</b>	GXEL00395
<b>Assay Type</b>	Quantitative
<b>Detection Method</b>	Colorimetric
<b>Sensitivity</b>	6.49 ng/mL
<b>Detection Range</b>	15.63 - 10, 000 ng/mL
<b>Recovery</b>	80-120%
<b>Sample Types</b>	Plasma, Serum
<b>Storage</b>	2-8°C. See Storage & Stability.
<b>Note</b>	For Research Use Only. Not for diagnostic procedures.

### Storage & Stability

When the kit was stored at the recommended temperature for 6 months, the signal intensity decreased by less than 20%. For unopened kits, if you want to prolong the storage time, please store the Standard, Detection A, Detection B and Microplate at -20 °C, the rest reagents should be store at 4°C. Store at 2-8°C on receipt; avoid repeated freeze-thaw cycles.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | [www.assaygenie.com](http://www.assaygenie.com) | [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com)  
 Datasheet: <https://store-h68l9z2lnx.mybigcommerce.com/content/Technical%20Manuals/ELISA/GXEL/GXEL00395.pdf> | Version 1.0, 09 July 2026