

Nanodisc Human OR2H1 Protein

HDFP009

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	HDFP009
Contents:	10µg Bradford Reagent: 1 vial (2ml)
Category:	Full Length Transmembrane Proteins
Synonyms:	6M1-16, dJ994E9.4, HS6M1-16, OLFR42A-9004-14, OLFR42A-9004.14/9026.2, OR2H6, OR2H8, OR6-2
Tag:	C-Flag Tag
Target:	OR2H1
Conjugation:	Unconjugated
Expression Host:	HEK293

Protein Data

Uniprot ID:	Q9GZK4
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Olfactory transduction
Formulation & Reconstitution:	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Molecular Weight:	The human full length OR2H1 Protein has a MW of 35.3 kDa
Storage & Shipping:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

-80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Store Bradford Reagent at Room Temperature for 1 year.

Usage: Research use only

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol