

Nanodisc Human CCR8 full length protein-MNP

HDFP130

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	HDFP130
Contents:	50µg Bradford Reagent: 1 vial (2ml)
Category:	Full Length Transmembrane Proteins
Synonyms:	CC-CKR-8, CCR-8, CDw198, CKRL1, CMKBR8, CMKBRL2, CY6, GPRCY6, TER1
Tag:	C-Flag Tag
Target:	CCR8
Conjugation:	Unconjugated
Expression Host:	HEK293

Protein Data

Uniprot ID:	P51685
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Chemokine signaling pathway, cytokine-cytokine receptor interaction
Formulation & Reconstitution:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Molecular Weight:	The human full length CCR8 Protein has a MW of 40.7 kDa
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at

-80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Store Bradford Reagent at Room Temperature for 1 year.

Usage: Research use only

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol