

## Nanodisc Human CLDN6 full length protein -MNP

HDFP132

### Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	HDFP132
<b>Contents:</b>	50µg Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	Full Length Transmembrane Proteins
<b>Synonyms:</b>	Claudin 6, Claudin-6, Skullin, Claudin6
<b>Tag:</b>	C-Flag Tag
<b>Target:</b>	CLDN6
<b>Conjugation:</b>	Unconjugated
<b>Expression Host:</b>	HEK293

### Protein Data

---

<b>Uniprot ID:</b>	P56747
<b>Protein Families:</b>	Transmembrane
<b>Protein Pathways:</b>	Cell adhesion molecules (CAMs), Leukocyte transendothelial migration, Tight junction
<b>Formulation &amp; Reconstitution:</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Molecular Weight:</b>	The human full length CLDN6 Protein has a MW of 23 kDa
<b>Storage &amp; Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

Store Bradford Reagent at Room Temperature for 1 year.

**Usage:** Research use only

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol