

## HDFP160

---

### Product Information

<b>Product SKU:</b>	HDFP160	<b>Expression Host:</b>	HEK293	<b>Size:</b>	10µg
<b>Target:</b>	5HT4R	<b>Tag:</b>	C-Flag Tag		

---

### Additional Information

<b>Conjugate:</b>	Unconjugated	<b>Uniprot ID:</b>	Q13639
<b>Molecular Weight:</b>	The human full length 5HT4R protein has a MW of 43.8kDa		

---

### Protein Information

**Background:** This gene is a member of the family of serotonin receptors, which are G protein coupled receptors that stimulate cAMP production in response to serotonin (5-hydroxytryptamine). The gene product is a glycosylated transmembrane protein that functions in both the peripheral and central nervous system to modulate the release of various neurotransmitters. Multiple transcript variants encoding proteins with distinct C-terminal sequences have been described. [provided by RefSeq, May 2010]

**Synonyms:** 5-HT4, 5-HT4R

**Protein Description:** Human 5HT4R full length protein-synthetic nanodisc

**Formulation:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with a pH below 6.5 or those containing high concentrations of divalent metal ions (greater than 5 mM) in subsequent experiments.

**Protein Pathways:** GPCRDB Class A Rhodopsin-like, Monoamine GPCRs.

**Protein Families:** GPCR, Transmembrane, Druggable Genome.

**Usage:** Research use only

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

**Contact Details | Dublin, Ireland**

**Email:** [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | **Web:** [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.