

## Bovine IL-17A Recombinant Protein

### Description

---

This high-purity Bovine IL-17A Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Protein Information

---

<b>SKU:</b>	KPRT0056
<b>Contents:</b>	Protein: 1 vial (5µg / 25µg / 100µg / 500µg(5X100µg)) Bradford Reagent: 1 vial (2ml)
<b>Source:</b>	Yeast
<b>Application:</b>	Cell Culture, ELISA, ELISpot, WB
<b>Molecular Weight:</b>	17.4 kDa
<b>Purity:</b>	>98% as visualized by SDS-PAGE analysis
<b>Amino Acid Sequence:</b>	AAIQSQDYTF RDISQKSLVL SGSYELRALH LNGQNMNQV VFRMSFVHGE ENSKKIPVVL CIKKNLYLS CVMKDGKPTL QLEMLDPKVY PKKKMEKRFV FNKTEIKGNV EFESSQFPNW YISTSQAEEM PVFLGNTKGG QDITDFIMES AS
<b>Homology:</b>	Acinonyx jubatus (cheetah), Felis catus (domestic cat), Leopardus geoffroyi (Geoffroy's cat), Lynx canadensis (Canada lynx), Lynx pardinus (Spanish lynx), Lynx rufus (bobcat), Panthera uncia (snow leopard), Prionailurus bengalensis (leopard cat), Prionailurus viverrinus (fishing cat)

## Preparation & Storage

---

**Stability & Storage:** Store Protein at -20°C for 1 Year.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.