

Equine CD40 Ligand (CD40LG) Recombinant Protein

Description

This high-purity Equine CD40 Ligand (CD40LG) Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

| | |
|-----------------------------|---|
| SKU: | KPRT0267 |
| Contents: | Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml) |
| Source: | Yeast |
| Application: | Cell Culture, ELISA, ELISpot, WB |
| Molecular Weight: | 16.0 kDa |
| Purity: | >98% as visualized by SDS-PAGE analysis |
| Amino Acid Sequence: | MQKGDQEPQI AAHVISEASS KTASVLQWAQ KGYTISNNL VTLENGKQLA VKRQGLYYIY AQVTFCSNRE ASGQAPFIAS LCLRSVSGSE RILLRAANTH SSSKPCGQQS IHLGGVFELQ PGASVFNVT DPSQVSHGTG FTSFGLLKL |
| Homology: | Equus asinus (ass), Equus caballus (horse), Equus przewalskii (Przewalski's horse) |

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.