

## Dolphin CXCL10 Recombinant Protein

### Description

---

This high-purity Dolphin CXCL10 Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Protein Information

---

|                             |   |
|-----------------------------|---|
| <b>SKU:</b>                 | KPRT0377  |
| <b>Contents:</b>            | Protein: 1 vial (5µg)<br>Bradford Reagent: 1 vial (2ml)   |
| <b>Source:</b>              | Yeast   |
| <b>Application:</b>         | Cell Culture, ELISA, ELISpot, WB  |
| <b>Molecular Weight:</b>    | 9.3 kDa   |
| <b>Purity:</b>              | >98% as visualized by SDS-PAGE analysis   |
| <b>Amino Acid Sequence:</b> | IPLSRTRCT CIKISDRPVN PRSLEKLEVI PASQSCPRVE IIATMKKNGE<br>KRCLNPESKT IKNLLKAISK ERSKRSQTQK EA  |
| <b>Homology:</b>            | Delphinapterus leucas (beluga whale), Globicephala melas (long-finned pilot whale), Lagenorhynchus obliquidens (Pacific white-sided dolphin), Neophocaena asiaeorientalis asiaeorientalis (Yangtze finless porpoise), Orcinus orca (killer whale), Phocoena sinus (vaquita), Tursiops truncatus (common bottlenose dolphin) |

## Preparation & Storage

---

**Stability & Storage:** Store Protein at -20°C for 1 Year.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.