

Multi-species CXCL12 (Ferret, Rabbit) Recombinant Protein

Description

This high-purity Multi-species CXCL12 (Ferret, Rabbit) Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT0387
Contents:	Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	8.6 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	GKPVLSYRC PCRFFESHVA RANVKHLKIL NTPNCALQIV ARLKSNNRQV CIDPKLKWIQ EYLEKALNKR FKM
Homology:	Dipodomys ordii (Ord's kangaroo rat), Dipodomys spectabilis (banner-tailed kangaroo rat), Heterocephalus glaber (naked mole-rat), Mustela erminea (ermine), Mustela putorius furo (domestic ferret), Ochotona princeps (American pika), Octodon degus (degu), Oryctolagus cuniculus (rabbit), Perognathus longimembris pacificus (Pacific pocket mouse), Sciurus carolinensis (gray squirrel)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.