

## Human/Monkey Midkine Recombinant Protein

### Description

---

This high-purity Human/Monkey Midkine Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Protein Information

---

<b>SKU:</b>	KPRT0521
<b>Contents:</b>	Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml)
<b>Source:</b>	Yeast
<b>Application:</b>	Cell Culture, ELISA, ELISpot, WB
<b>Molecular Weight:</b>	13.4 kDa
<b>Purity:</b>	>98% as visualized by SDS-PAGE analysis
<b>Amino Acid Sequence:</b>	VAKKKDKVKK GPGSECAEW AWGPCTPSSK DCGVGFREGT CGAQTQRIRC RVPCNWKKEF GADCKYKFEN WGACDGGTGT KVRQGTLLKKA RYNAQCQETI RVTKPCTPKT KAKAKAKKGK GKD
<b>Homology:</b>	Chlorocebus sabaeus (green monkey), Gorilla gorilla gorilla (western lowland gorilla), Homo sapiens (human), Macaca mulatta (Rhesus monkey), Mandrillus leucophaeus (drill), Papio anubis (olive baboon), Pongo abelii (Sumatran orangutan), Rhinopithecus bieti (black snub-nosed monkey), Rhinopithecus roxellana (golden snub-nosed monkey), Theropithecus gelada (gelada)

## Preparation & Storage

---

**Stability & Storage:** Store Protein at -20°C for 1 Year.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.