

## Bovine IL-25 Recombinant Protein

### Description

---

This high-purity Bovine IL-25 Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Protein Information

---

<b>SKU:</b>	KPRT0554			
<b>Contents:</b>	Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml)			
<b>Source:</b>	Yeast			
<b>Application:</b>	Cell Culture, ELISA, ELISpot, WB			
<b>Molecular Weight:</b>	17.5 kDa			
<b>Purity:</b>	>98% as visualized by SDS-PAGE analysis			
<b>Amino Acid Sequence:</b>	LCLGSQRCCT	HWPSCCPSKG	QNSTEEWLKW	NSVLMPPPET
	TSLAHHPEPC	SSSRDGPLNS	RSIAPWRYEL	DRDLNRLPQD
	LYHARCLCPH	CVSLQTGSHM	DPLGNSELY	HNQTVFYRRP
	CPGQQGDHGG	YCLERRLYRV	SLACVCVRPR	VMA
<b>Homology:</b>	Bison bison bison (American buffalo), Bos indicus (zebu), Bos mutus (wild yak), Bos taurus (cattle)			

## Preparation & Storage

---

**Stability & Storage:** Store Protein at -20°C for 1 Year.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.