

Chicken and Turkey TGF-alpha Recombinant Protein

Description

This high-purity Chicken and Turkey TGF-alpha Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

| | | | | |
|-----------------------------|---|------------|------------|--------------------------|
| SKU: | KPRT0555 | | | |
| Contents: | Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml) | | | |
| Source: | Yeast | | | |
| Application: | Cell Culture, ELISA, ELISpot, WB | | | |
| Molecular Weight: | 5.7 kDa | | | |
| Purity: | >98% as visualized by SDS-PAGE analysis | | | |
| Amino Acid Sequence: | VVSHFNKCPD | SHTQYCFHGT | CRFLVQEEKP | ACVCHSGYVG VRCEHADLLA |
| Homology: | Coturnix japonica (Japanese quail), Gallus gallus (chicken), Meleagris gallopavo (turkey), Numida meleagris (helmeted guineafowl) | | | |

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.