

Porcine IL-25 Recombinant Protein

Description

This high-purity Porcine IL-25 Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT0574		
Contents:	Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml)		
Source:	Yeast		
Application:	Cell Culture, ELISA, ELISpot, WB		
Molecular Weight:	17.1 kDa		
Purity:	>98% as visualized by SDS-PAGE analysis		
Amino Acid Sequence:	WSQKGCPHWP	NCCPSKGQNP	IEEWLKQNAV
	TPHPESCKAS	EDGPLNSRSI	APWRYELDRD
	ARCLCPHCVS	LRTGSHMDPL	GNSELLYHNQ
	QQGAHDSYCL	EPRLYQVSLA	CVCVRPRVMA
Homology:	Sus scrofa (pig)		

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.