

Chicken Annexin V Recombinant Protein

Description

This high-purity Chicken Annexin V Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT0576
Contents:	Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	36.2 kda
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	MAKYTRGTVT AFSPFDARAD AEALRKAMKG MGTDEETILK ILTSRNNAQR QEIASAFKTL FGRDLVDDLK SELTGKFETL MVSLMRPARI FDAHALKHAI KGAGTNEKVL TEILASRTPA EVQNIKQVYM QEYEANLEDK ITGETSGHFQ RLLVLLQAN RDPDGRVDEA LVEKDAQVLF RAGELKWGTD EETFITILGT RSVSHLRRVF DKYMTISGFQ IEETIDRETS GDLEKLLAV VKCIRSVPAY FAETLYYSMK GAGTDDDTLI RVMVSRSEID LLDIRHEFRK NFAKSLYQMI QKDTSGDYRK ALLLLCGGDD E
Homology:	Gallus gallus (chicken)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.