

Bovine CD27 Recombinant Protein

Description

This high-purity Bovine CD27 Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT0592			
Contents:	Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml)			
Source:	Yeast			
Application:	Cell Culture, ELISA, ELISpot, WB			
Molecular Weight:	17.7 kDa			
Purity:	>98% as visualized by SDS-PAGE analysis			
Amino Acid Sequence:	TPGPKSCPEK	HYWAQGGWCC	QMCEPGTFLV	KDCEQHREAA
	QCNPCTPGVS	FMPDHHSRPH	CESCRHCNSG	LLIRNCTLTA
	NSECACPEGQ	QCRDKDCMEC	DGPAQAPGPH	PQPSHLPYAE
	EIAESRTDRH	TQTLANSRWL	PAPTLSTHWS	PQRSLCSANC
Homology:	Bos indicus (Zebu), Bos taurus (cattle)			

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.