

Equine Annexin V Recombinant Protein

Description

This high-purity Equine Annexin V Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT0632
Contents:	Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	36.0 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	MAQVLKGTVT DFPGFDERAD AEVLRKAMKG LGTDEESILT LLTSRSNAQR QEIAVAFKTL FGRDLLDDLK SELTGKFEKL IVALMKPSRL YDAYELKHAL KGAGTNEKVL TEIIASRTPA ELRAIKQVYE EEYGSNLEDD VVADTSGFYQ RMLVVLLQAN RDPDAGIDEA QVEQDAQTLF QAGELKWGTD EEKFITIFGT RSVSHLRRVF DKYMTISGFQ IEETIDRETS GNLEQLLAV VKSIRSIPAY LAETLYYAMK GAGTDDHTLI RVVVSRSSEID LFNIRKEFRK NFATSLYSMI KSDTSGDYKK ALLLLCGGED D
Homology:	Equus caballus (horse), Equus przewalskii (Przewalski's horse)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.