

Bovine GM-CSF Recombinant Protein

Description

This high-purity Bovine GM-CSF Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT0871
Contents:	Protein: 1 vial (25µg, 100µg, 500µg(5X100µg)) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	14.3 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	APTRPPNTAT RPWQHVDAIK EALSLLNHSS DTDVMNDTE VVSEKFDSQE PTCLQTRLKL YKNGLQGLT SLMGSLTMMMA THYEKHCPT PETSCGTQFI SFKNFKEDLK EFLFIIPFDC WEPAQK
Homology:	Bos taurus (cattle)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.