

Canine Annexin V Fluorescein Recombinant Protein

Description

This high-purity Canine Annexin V Fluorescein Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT0964
Contents:	Protein: 1 vial (10µg(100tests)) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	MAQVLKGTVT PFPGFDERAD AEVLRKAMKG LGTDEESILT LLTSRSNAQR QEIAAAFCTL FGRDLLDDLK SELTGKFEKL IVALMKPSRL YDAYELKHAL KGAGTNEKVL TEIIASRTPE ELRAIKQVYE EEYGSSLEDD VVGDTSGYYQ RMLVVLLQAN RDPDAGIDEA QVEQDAQALF QAGELKWGTD EEKFITIFGT RSVSHLRRVF DKYMTISGFQ IEETIDRETS GNLEQLLAV VKSIRSIPAY LAETLYYAMK GAGTDDHTLI RVVVSRSSEID LFNIRKEFRK NFATSLYSMI KGDTSGDYKK ALLMLCGGED
Homology:	Canis lupus familiaris (dog)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.