

Ovine IL-4 Recombinant Protein

Description

This high-purity Ovine IL-4 Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT1130
Contents:	Protein: 1 vial (25µg, 100µg, 500µg(5X100µg)) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	12.6 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	HKCDITL E E I I K T L N I L T S R K N S C M E L P V A D V F A A P K N A T E K E T F C R A G I E L R R I Y R S H M C L N K F L G G L D R N L S S L A S K T C S V N E A K T S T S T L R D L L E R L K T I M K E K Y S K C
Homology:	Budorcas taxicolor (takin), Ovis aries (sheep)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.