

Chicken LIF Recombinant Protein

Description

This high-purity Chicken LIF Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT1395
Contents:	Protein: 1 vial (25µg, 100µg, 500µg(5X100µg)) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	21.0 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	RALLGTSSAC PTNGLCRANV LEQTRRQVAL LNATAQDLFS LYLKCQGEPP SSESDRLCSP SGIFFPPFHV NRITTERKEVM VAMYKLFAFL NASLGNITRD QEELNPMAKE LLNRLHNNTK TTRGLISNLT CLLCKHYNIF QVDVSYGESS KDKSAFKKKQ QGCQVLRKYV QVIAQAARVL LPHLSPA
Homology:	Gallus gallus (chicken)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.