

Chicken PDGF-BB Recombinant Protein

Description

This high-purity Chicken PDGF-BB Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

| | |
|-----------------------------|---|
| SKU: | KPRT1397 |
| Contents: | Protein: 1 vial (25µg, 100µg, 500µg(5X100µg)) Bradford Reagent: 1 vial (2ml) |
| Source: | Yeast |
| Application: | Cell Culture, ELISA, ELISpot, WB |
| Molecular Weight: | 12.4 kDa (monomer) |
| Purity: | >98% as visualized by SDS-PAGE analysis |
| Amino Acid Sequence: | SLDALAAEP AVLAECKTRT VVFEISRDMV DSTNANFVWW PPCVEVQRCS GCCNNRNVQC RPMQIRVRHV QVNKIEFFQR KPIFKKVIVP LEDHVQCRCE VSRPPPRS |
| Homology: | Gallus gallus (chicken) |

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.