

Cynomolgus Monkey M-CSF Recombinant Protein

Description

This high-purity Cynomolgus Monkey M-CSF Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT1612
Contents:	Protein: 1 vial (25µg, 100µg, 500µg(5X100µg)) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	18.4 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	EEVSEYCSHM IGSGHLQSLQ RLIDSQMETS CQITFEFVDQ EQLKDPVCYL KKAFLLVQDI MEDTMRFRDN TPNAIAIVQL QELSLRLKSC FTKDYEEHDK ACVRTFYETP LQLEKVKNV FNETKNLLDK DWNIFSKNCN NSFAECSGQG HERQREGS
Homology:	Cercocebus atys (sooty mangabey), Chlorocebus sabaeus (green monkey), Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey), Papio anubis (olive baboon), Ptilocolobus tephrosceles (Ugandan red Colobus), Rhinopithecus roxellana (golden snub-nosed monkey), Theropithecus gelada (gelada)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.