

## Human TNF-R1 Recombinant Protein

### Description

---

This high-purity Human TNF-R1 Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Protein Information

---

<b>SKU:</b>	KPRT1807			
<b>Contents:</b>	Protein: 1 vial (5µg, 25µg) Bradford Reagent: 1 vial (2ml)			
<b>Source:</b>	Yeast			
<b>Application:</b>	Cell Culture, ELISA, ELISpot, WB			
<b>Molecular Weight:</b>	20.5 kDa			
<b>Purity:</b>	>98% as visualized by SDS-PAGE analysis			
<b>Amino Acid Sequence:</b>	LVPHLGDREK	RDSVCPQGKY	IHPQNNSICC	TKCHKGTYLY
	NDCPGPGQDT	DCRECESGSF	TASENHLRHC	LSCSKCRKEM
	GQVEISSCTV	DRDTVCGCRK	NQYRHYWSEN	LFQCFNCSLC
	LNGTVHLSCQ	EKQNTVCTCH	AGFFLRENEC	VSCSNCKKSL
	ENVKGTEDSG	TT		ECTKLCLPQI
<b>Homology:</b>	Gorilla gorilla gorilla (western lowland gorilla), Homo sapiens (human), Pongo abelii (Sumatran orangutan)			

## Preparation & Storage

---

**Stability & Storage:** Store Protein at -20°C for 1 Year.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.