

Bovine TNF alpha Biotinylated Recombinant Protein

Description

This high-purity Bovine TNF alpha Biotinylated Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT1817
Contents:	Protein: 1 vial (10µg) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	LRSSSQASSN KPAHVADI NSPGQLRWWD SYANALMANG VKLEDNQLVV PADGLYLIYS QVLFGRGQCP STPLFLTHI SRIAVSYQTK VNILSAIKSP CHRETPWAE AKPWYEPIYQ GGVFQLEKGD RLSAEINLPD YLDYAESGQV YFGIIAL
Homology:	Bison bison bison (American buffalo), Bos indicus (zebu), Bos javanicus (banteng), Bos mutus (wild yak), Bos taurus (cattle), Bubalus bubalis (water buffalo), Bubalus kerabau (carabao)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.