

Bovine VEGF-A Biotinylated Recombinant Protein

Description

This high-purity Bovine VEGF-A Biotinylated Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT1827
Contents:	Protein: 1 vial (10µg) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	APMAEGGQKP HEVVKFMDVY QRSFCRPIET LVDIFQEYPD EIFIFKPSC VPLMRCGGCC NDESLECVPT EEFNITMQIM RIKPHQSQHI GEMSFLQH NK CECRPKKDKA RQENPCGPCS ERRKHLFVQD PQTCKCCKN TDSRCKARQL ELNERTCRCD KPRR
Homology:	Bison bison bison (American buffalo), Bos taurus (cattle), Bubalus bubalis (water buffalo), Bubalus kerabau (carabao)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.