

## Cynomolgus Monkey CRP Recombinant Protein

### Description

---

This high-purity Cynomolgus Monkey CRP Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Protein Information

---

|                             |  |
|-----------------------------|--|
| <b>SKU:</b>                 | KPRT2151   |
| <b>Contents:</b>            | Protein: 1 vial (5µg, 25µg)<br>Bradford Reagent: 1 vial (2ml)  |
| <b>Source:</b>              | Yeast  |
| <b>Application:</b>         | Cell Culture, ELISA, ELISpot, WB   |
| <b>Molecular Weight:</b>    | 23.2 kDa   |
| <b>Purity:</b>              | >98% as visualized by SDS-PAGE analysis  |
| <b>Amino Acid Sequence:</b> | QTDMSMKAFV FPKESDNSYV TLKARLTKPL KAFTVCLHFY TELSSTRGYS<br>IFSYATKRQN NEILIFWSKD IGYSFTVGGG EILFEVPEVT VAPVHICTSW<br>ESASGIVEFW VDGKPRARKS LKRGYTVGED ASIILGQEQD SFGGSFETQQ<br>SLVGDIGNVN MWDFVLSPDE ISTVYRGGTF SPSVLYWRAL KYEVQGEVFI<br>KPQLWS |
| <b>Homology:</b>            | Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey)  |

## Preparation & Storage

---

**Stability & Storage:** Store Protein at -20°C for 1 Year.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.