

## Human IL-11 Recombinant Protein

### Description

---

This high-purity Human IL-11 Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Protein Information

---

<b>SKU:</b>	KPRT2153
<b>Contents:</b>	Protein: 1 vial (5µg, 25µg) Bradford Reagent: 1 vial (2ml)
<b>Source:</b>	Yeast
<b>Application:</b>	Cell Culture, ELISA, ELISpot, WB
<b>Molecular Weight:</b>	18.4 kDa
<b>Purity:</b>	>98% as visualized by SDS-PAGE analysis
<b>Amino Acid Sequence:</b>	RVSPDPRAEL DSTVLLTRSL LADTRQLAAQ LRDKFPADGD HNLDLPTLA MSAGALGALQ LPGVLTRLRA DLLSYLRHVQ WLRRAGGSSL KTLEPELGT LARLDRLRR LQLLMSRLAL PQPPDPPAP PLAPPSSAWG GIRAAHAILG GLHLTDWAV RGLLLLKTRL
<b>Homology:</b>	Gorilla gorilla gorilla (western lowland gorilla), Homo sapiens (human), Hylobates moloch (silvery gibbon), Nomascus leucogenys (northern white-cheeked gibbon), Pan paniscus (pygmy chimpanzee), Pan troglodytes (chimpanzee), Pongo abelii (Sumatran orangutan), Pongo pygmaeus (Bornean orangutan)

## Preparation & Storage

---

**Stability & Storage:** Store Protein at -20°C for 1 Year.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.