

Porcine IL-12p70 Recombinant Protein

Description

This high-purity Porcine IL-12p70 Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT2155
Contents:	Protein: 1 vial (5µg) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	56.5 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	IWELEKNVYV VELDWYPNAP GEMVVLTCNT PEEDGITWTS DQSSEVLGTG KTLTIHVKEF GDAGQYTCRK GGAVLSQSLI LLHKKEDGIW STDILKDQKE PKNKSFLKCE AKNYSGRFTC WWLTAISTDL KFSVKSSRGS TDPRGVTCGT ATLSEDLGEY KKYRVECQEG SACPAAEESL PIEVVLEAVH KLKYENYTSS FFIRDIKPD PPKNLQLNPL KNSRHVEISW EYPDTWSTPH SYFSLMFGVQ VQGKKNKREK DKLFTDQISA KVTCHKDANI RVQARDRYYS SSWSEWASVS CN (302) LINKER RSLPATTAGP GMFKCLNHSQ NLLKAVSNTL QKAKQTLEFY SCTSEEIDHE DITKDKTSTV EACLPLELAT NESCLAARET SLITNGNCLT SGKTSFMTTL CLSSIEDLK MYHVEFQAMN AKLLMDPKRQ IFLDQNMLTA ITELQALNF NSETVPQKPS LEELDFYKTK IKLCILLHAF RIRAVTIDRM MSYLNSS
Homology:	Sus scrofa (pig)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.