

Feline PD-L1 Recombinant Protein

Description

This high-purity Feline PD-L1 Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT2204
Contents:	Protein: 1 vial (5µg, 25µg, 100µg, 500µg(5X100µg)) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	24.6 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	FTITVSKDLY VVEYGSNVTM ECRFPVEEQL DLVSLIVYWE MEDKKIIQFV QGKEDLKVQH RSYSQRAQLL KDQLFLGKAA LQITNVTLED AGVYCCLIGY GGADYKRITL KVHAPYRKIN QRISVDPVTS EHELMCQAEG YPTAEVIWTN SAHQVLNGKT IISVSNMETK LFNVTSTLRI NTTANEIFYC TFLQRSSPEG NSTAELVIPE PFLVPANER
Homology:	Acinonyx jubatus (cheetah), Felis catus (domestic cat), Lynx canadensis (Canada lynx), Lynx rufus (bobcat), Puma concolor (puma)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.