

Ferret/Mink PDGF-BB Recombinant Protein

Description

This high-purity Ferret/Mink PDGF-BB Recombinant Protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU:	KPRT2227
Contents:	Protein: 1 vial (5µg, 25µg) Bradford Reagent: 1 vial (2ml)
Source:	Yeast
Application:	Cell Culture, ELISA, ELISpot, WB
Molecular Weight:	18.3 kDa
Purity:	>98% as visualized by SDS-PAGE analysis
Amino Acid Sequence:	SLGSPTVAEP AMIAECKTRT EVFEISRRLI DRTNANFLVW PPCVEVQRCS GCCNNRNVQC RPTQVQLRPV QVRKIEIVRK KPTFKKATVT LEDHLACKCE TVVATRPVTR SPGSTQEQRA AKTPQTRVTI RTVRVRRPPK GKHRKFKHHTH DKKALKETLG A
Homology:	Enhydra lutris kenyoni (northern sea otter), Lontra canadensis (Northern American river otter), Lutra lutra (Eurasian river otter), Mustela lutreola (European mink), Mustela nigripes (black-footed ferret), Mustela putorius furo (domestic ferret), Neogale vison (American mink)

Preparation & Storage

Stability & Storage: Store Protein at -20°C for 1 Year.
Store Bradford Reagent at Room Temperature for 1 Year.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.