

RPCB0075

Product Information

Product SKU:	RPCB0075	Gene ID:	131450	Size:	100µg
Tag:	C-His&Avi	Reactivity:	Human		

Additional Information

Expression Host:	-	Swissprot:	AAQ19772
Purity:	> 95% by Tris-Bis PAGE;> 95% by SEC-HPLC		

Protein Information

Background: CD200Fc, a chimeric molecule including the extracellular domain of CD200 and a murine IgG2a Fc region, regulates immune responses following engagement of a cell surface receptor, CD200R, expressed on cells of the myeloid and T cell lineage. A recent report focused attention on a family of CD200Rs, but concluded that only one member used CD200 as its ligand.

Protein Description: High quality, high purity and low endotoxin recombinant Biotinylated Recombinant Human CD200R1 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80°C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.