

Recombinant Human NKAT-2/KIR2DL3/CD158b2 Protein

RPCB0389



Product Information

Product SKU: RPCB0389 **Gene ID:** 3804 **Size:** 20µg
Tag: C-hFc&His **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** AAB36590.1
Purity: > 95% by SDS-PAGE.

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human NKAT-2/KIR2DL3/CD158b2 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 0.1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.