

RPCB0827

Product Information

Product SKU:	RPCB0827	Gene ID:	3845	Size:	10µg
Tag:	N-His	Reactivity:	Human		

Additional Information

Expression Host:	E. coli	Swissprot:	P01116
Purity:	> 85% by SDS-PAGE.		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human GTPase KRas/KRAS Protein , tested reactivity in E. coli and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 1.0 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4. Contact us for customized product form or formulation.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.