

Recombinant Human Serpin D1/Heparin cofactor 2 Protein

RPCB0892

Product Information

Product SKU: RPCB0892 **Gene ID:** 3053 **Size:** 20µg
Tag: C-His **Reactivity:** Human

Additional Information

Expression Host: HEK293 cells **Swissprot:** P05546
Purity: > 90% by SDS-PAGE.

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Serpin D1/Heparin cofactor 2 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 0.1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of 50mM Tris,150mM NaCl, pH7.5.Contact us for customized product form or formulation.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.