

RPCB1001

Product Information

Product SKU:	RPCB1001	Gene ID:	54790	Size:	20µg
Tag:	N-his	Reactivity:	Human		

Additional Information

Expression Host:	-	Swissprot:	Q6N021-1
Purity:	-		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human TET2 Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method.

Formulation: Supplied as a 0.22 µm filtered solution in 50mM HEPES,300mM NaCl,1mM TCEP,10% Glycerol,pH 8.0.

Storage: Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.