

# Recombinant Human Calgranulin-A&Calgranulin-B/S100A8&S100A9 Protein

RPCB1345



## Product Information

<b>Product SKU:</b>	RPCB1345	<b>Gene ID:</b>	6279&6280	<b>Size:</b>	100µg
<b>Tag:</b>	C-His	<b>Reactivity:</b>	Human		

## Additional Information

<b>Expression Host:</b>	E.coli	<b>Swissprot:</b>	P05109&P06702
<b>Purity:</b>	> 95% by Tris-Bis PAGE;> 95% by SEC-HPLC		

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant Human Calgranulin-A&Calgranulin-B/S100A8&S100A9 Protein , tested reactivity in E.coli and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 1 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of 2mM DTT, PBS (pH 7.2).

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.