

RPCB1354

Product Information

Product SKU:	RPCB1354	Gene ID:	6772	Size:	10µg
Tag:	C-His	Reactivity:	Human		

Additional Information

Expression Host:	Baculovirus-Infected Sf9 Cells	Swissprot:	P42224-2
Purity:	> 92% by SDS-PAGE.		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human STAT1 Protein , tested reactivity in Baculovirus-Infected Sf9 Cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 1.0 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.