

RPCB1447

Product Information

Product SKU:	RPCB1447	Gene ID:	928938	Size:	10µg
Tag:	C-His	Reactivity:	Monkeypox virus		

Additional Information

Expression Host:	E. coli	Swissprot:	Q8V4Z7
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Monkeypox virus L1R Protein , tested reactivity in E. coli and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: <0.1 EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.