

RPCB1478

Product Information

Product SKU:	RPCB1478	Gene ID:	43740568	Size:	100µg
Tag:	C-His&Flag	Reactivity:	SARS-COV-2		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	QHD43416.1
Purity:	> 95% by Tris-Bis PAGE;> 95% by SEC-HPLC		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant SARS-COV-2 Spike S1 NTD Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 1 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.